

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

CRS (CysteinyI tRNA synthetase)

产品货号: TP01720

第三版

别名:

描述: CRS is a class 1 aminoacyl-tRNA synthetase, cysteinyI-tRNA synthetase. Each of the twenty aminoacyl-tRNA synthetase catalyzes the aminoacylation of a specific tRNA or tRNA isoaccepting family with the cognate amino acid.

配方: 1X PBS (pH7.2) 200mM NaCl 20% Glycerol

分子量: 85kDa

序列:

纯度: > 95% by HPLC

浓度: 0.50mg/ml

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.