



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

CTHRC1, 31-243 aa, Human, His tag, E.coli

产品货号: TP01767

第三版

别名: Collagen triple helix repeat-containing protein 1 isoform 1, NMTC1

描述: CTHRC1 also known as Collagen triple helix repeat-containing protein 1. This protein may act as a negative regulator of collagen matrix deposition. Recombinant human CTHRC1 protein, was expressed in E.coli

配方: Liquid. In 20mM Tris-HCl buffer (pH 8.0) containing 10% glycerol

分子量: 25.3kDa (234aa)

序列:

MGSSHHHHHHSSGLVPRGSHMSEIPKGKQKAQLRQREVVDLYNGMCLQGPAVPGRDGSPGANGIPGTPGIPGRDGFK  
GEKGECLRESFEESWTPNYKQCSWSSLNYGIDLGKIAECTFTKMRSNSALRVLFSGSLRLKCRNACCQRWYFTFNGAECESG  
PLPIEAIYLDQGSPEMNSTINIHRTSSVEGLCEGIGAGLVDVAIWVGTCSDYKGDASTGWNSVSRIIIIEELPK

纯度: > 95% by HPLC

浓度: 1mg/ml (determined by Bradford assay)

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.