

本公司提供的电版本说明书仅供参考，实验请以收到的纸质手册为准。

CTNNBIP1, 1-81aa, Human, His tag, E.coli

产品货号: TP01772

第三版

别名: Beta-catenin-interacting protein 1, ICAT

描述: CTNNBIP1 binds CTNNB1 and prevents interaction between CTNNB1 and TCF (T-cell transcription factor) family members. This protein is a negative regulator of the Wnt signaling pathway.

配方: Liquid. 20mM Tris-HCl buffer (pH8.0) containing 20% glycerol, 2mM DTT, 100mM NaCl

分子量: 11.3kDa (101a) confirmed by MALDI-TOF

序列:

MGSSHHHHHHSSGLVPRGSHMNREGAPGKSPEEMYIQQKVRVLLMLRKMGSNLTASEEEFLRTYAGVVNSQLSQLPPHS
IDQGAEDVVMAFSRSETEDRRQ

纯度: > 95% by HPLC

浓度: 1 mg/ml (determined by Bradford assay)

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.