

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

CXCL2, 35-107 aa, Human, His-tagged, Recombinant, E.coli

产品货号: TP01798

第三版

别名:C-X-C motif chemokines 2, MIP2-alpha, Gro-beta, CINC-2a, GRO2, GROb, MIP2, MIP2A, SCYB2

描述: CXCL2 is a small cytokine belonging to the CXC chemokine family that is also known as Chemokine (C-X-C motif) ligand 2, and macrophage inflammatory protein 2-alpha (MIP2-alpha). CXCL2 is secreted by activated monocytes, neutrophils and macrophages and expressed at sites of inflammation.

配方: Liquid. In 20 mM Tris-HCl buffer (pH 8.0)

分子量: 10.1 kDa (94aa)

序列:

MGSSHHHHHHSSGLVPRGSHMAPLATELRCQCLQTLQGIHLKNIQSVKVKSPGPHCAQTEVIATLKNQKACLNPASP  
MVKKIIEKMLKNGKSN

纯度: > 95% by HPLC

浓度: 1.0 mg/ml (determined by Bradford assay)

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.