

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

DAAO, 1-347 aa, Human, His-tagged, Recombinant, E.coli

产品货号: TP01828

第三版

别名:D-amino acid oxidase, DAMOX, DAO, OXDA

描述:D-amino-acid oxidase (DAAO) is a peroxisomal enzyme which uses flavin adenine dinucleotide (FAD) as a cofactor and oxidizes D-amino acids to the corresponding imino acids, producing ammonia and hydrogen peroxide. Its substrates include a wide variety of D-amino acids, but it is inactive on the naturally occurring L-amino acids.

配方:Liquid. In 20 mM Tris-HCl buffer (pH 8.0) containing 20% glycerol 1mM DTT

分子量:41.6 kDa (367aa)

序列:

```
MGSSHHHHHHSSGLVPRGSHMRVVVIGAGVIGLSTALCIHERYHSQLPLDIKVYADRFTPLTTTDDVAAGLWQPYLSDP  
NNPQEADWSQQTDFDYLLSHVHSPNAENLGLFLISGYNLFHEAIPDPSWKDTVLGFRKLTPRELDMFPDYGYGWFHTSLILE  
GKNYLQWLTERLTERGVKFFQRKVESFEEVAREGADVIVNCTGVWAGALQRDPLLQPGRGQIMKVDAPWMKHFILTH  
DPERGIYNSPYIIPGTQTVTLGGIFQLGNWSELNNIQDHNTIWEGCCRLEPTLKNARIIGERTGFRPVRPQIRLEREQLRTGP  
SNTEVIHNYGHGGYGLTIHWGCALEAAKLFGRILEEKKLSRMPPSHL
```

纯度:> 95% by HPLC

浓度:0.5 mg/ml (determined by Bradford assay)

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.