

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

DAP, 1-102aa, Human, His tag, E.coli

产品货号: TP01829

第三版

别名: death-associated protein, DAP-1, DAP1

描述: DAP is a basic, proline-rich, 15-kD protein. The protein acts as a positive mediator of programmed cell death that is induced by interferon-gamma. Recombinant human DAP protein, fused to His-tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography techniques.

配方: Liquid. In 20mM Tris-HCl buffer (pH 8.0) containing 0.15M NaCl, 10% glycerol, 1mM DTT

分子量: 13.6kDa(125aa) confirmed by MALDI-TOF

序列:

MGSSHHHHHHSSGLVPRGSHMGSMSSPPEGKLETKAGHPPAVKAGGMRIVQKHPHTGDTKEEKDKDDQEWESPSPPKP  
TVFISGVIARGDKDFPPAAAQVAHQKPHASMDKHPSRPTQHIQQPRK

纯度: > 95% by HPLC

浓度: 0.25 mg/ml (determined by BRADFORD assay)

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.