

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

DBNDD1, 1-158aa, Human, His tag, E.coli

产品货号: TP01832

第三版

别名:Dysbindin domain-containing protein 1 isoform 1.

描述:DBNDD1 also known as dysbindin domain-containing protein 1, belongs to the dysbindin family. Recombinant human DBNDD1 protein, fused to His-tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography.

配方:Liquid. 20mM Tris-HCl buffer (pH8.0) containing 10% glycerol, 0.1M NaCl, 1mM DTT

分子量:19.6 kDa (182aa) confirmed by MALDI-TOF

序列:

MGSSHHHHHHSSGLVPRGSHMGSHMEPPEGAGTGEIVKEAEVPAALGVPAQGTGDNGHTPVVEEVGGIPVPAPGLLQ
VTERRQLSSVSSLEVHFDLLDLTELTDMSDQELAEVFADSDDENLNTESPAGLHPLPRAGYLRSWTRTRAEQSHEKQPL
GDPERQATVLDTFLTVERPQED

纯度:> 95% by HPLC

浓度:0.5 mg/ml

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C 或者-70 ° C. 避免反复冻融.