

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

DCTN2, 1-406aa, Human, His tag, E.coli

产品货号: TP01843

第三版

别名: Dynactin subunit 2 isoform 1, Dynactin 2 (p50), DCTN50, DYNAMITIN, RBP50

描述: DCTN2 is a 50-kD subunit of dynactin, a macromolecular complex consisting of 10-11 subunits ranging in size from 22 to 150 kD. Dynactin binds to both microtubules and cytoplasmic dynein. It is involved in a diverse array of cellular functions, including ER-to-Golgi transport, the centripetal movement of lysosomes and endosomes, spindle formation, chromosome movement, nuclear positioning, and axonogenesis. This subunit is present in 4-5 copies per dynactin molecule. It contains three short alpha-helical coiled-coil domains that may mediate association with self or other dynactin subunits. It may interact directly with the largest subunit (p150) of dynactin and may affix p150 in place. Recombinant human DCTN2 protein, fused to His-tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography techniques.

配方: Liquid. In 20mM Tris-HCl(pH8.0) containing 20% glycerol, 0.15M NaCl, 1mM DTT

分子量: 47.2kDa (429aa) confirmed by MALDI-TOF

序列:

MGSSHHHHHSSGLVPRGSHMGSMADPKYADLPGIARNEPDVYETSDLPEDDQAEFDFAQEELELTSTSV EHIIVNPNA
AYDKFKDKRVGKGLDFSDRIGKTKRTGYESGEYEMLGEGLGVKETPQQKYQRLLHEVQELTTEVEKIKTTVKESATEEKL
TPVLLAKQLAALKQQLVASHLEKLLGPDAAINLTPDGALAKRLLLQLEATKNSKGGSGGKTTGTPPDSSLVTYELHSRPE
QDKFSQAAKVAELEKRLTELETAVRCDQDAQNPLSAGLQGA CLMETVELLQAKVSALDLAVLDQVEARLQSVLGKVNEI
AKHKASVEDADTQSKVHQLYETIQRWSP IASTLPELVQRLVTIKQLHEQAMQFGQLLTHLDTTQQMIANSLKDNTLLT
QVQTTMRENLATVEGNFASIDERMK KLGK

纯度: > 95% by HPLC

浓度: 0.5 mg/ml (determined by Bradford assay)

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C 或者-70 ° C. 避免反复冻融.