

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

DCUN1D1, 1-259aa, Human, His tag, E.coli

产品货号: TP01846

第三版

别名:DCN1-like protein 1, DCUN1L1, RP42, SCCRO, SCRO, Tes3

描述:DCUN1D1, also known as DCN1-like protein 1, is involved in the malignant transformation of squamous cell lineage. This protein regulates Gli1, a key regulator of the hedgehog (HH) pathway that plays an important role in development, maintenance, and regeneration of almost all adult tissues.

配方:Liquid. 20mM Tris-HCl buffer (pH8.0) containing 10% glycerol 0.1 M NaCl

分子量:32.2 kDa (279aa), confirmed by MALDI-TOF

序列:

MGSSHHHHHHSSGLVPRGSHMNKLKSSQKDKVRQFMIFTQSSEKTAVSCLSQNDWKLDVATDNFFQNPELYIRESVKGS
LDRKKLEQLYNRYKDPQDENKIGIDGIQQFCDDLALDPASISVLIIAWKFRAATQCEFSKQEFMDGMTELGCDSEKLAQI
PKMEQELKEPGRFKDFYQFTFNFAKNPGQKGLDLEMAIAYWNLVLNGRFKFLDLWNKFLLEHHKRSIPKDTWNLLDFST
MIADDMSNYDEEGAWPVLIDDFVEFARPQIAGTKSTTV

纯度:> 95% by HPLC

浓度:1.0 mg/ml (determined by Bradford assay)

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.