

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

DCUN1D4, 1-292aa, Human, His tag, E.coli

产品货号: TP01847

第三版

别名:DCN1, defective in cullin neddylation 1, domain containing 4, DCN1-like protein 4 isoform 1

描述:DCUN1D4 contains 1 DCUN1 domain. Recombinant human DCUN1D4 protein, fused to His-tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography techniques.

配方:Liquid. In 20mM Tris-HCl buffer (pH 8.0) containing 0.2M NaCl, 50% glycerol, 2mM DTT, 2mM EDTA

分子量:36.5kDa (315aa) confirmed by MALDI-TOF

序列:

MGSSHHHHHHSSGLVPRGSHMGSMHSDAAAVNFQLNSHLSTLANIHKIYHTLNKLNLTEDIGQDDHQTGSLRSCSSSDC
FNKVMPPRKKRRPASGDDLAKKSRHDSMYRKYDSTRIKTEEEAFSSKRCLEWYFYAGTDDVVGPEGMEKFCEDIGVEP
ENVVMLVLAWKLDQNMGYFTLQEWLKGMTSLQCDTTEKLRNTLDYLRSLNDSTNFKLIYRYAFDFAREKDQQRSLDIN
TAKCMLGLLLGIWPLFPVFHQFLEQSKYKVINKDQWCNVLEFSRTINLDLSNYDEDGAWPVLLDEFVEWYKDKQMS

纯度:> 95% by HPLC

浓度:0.25 mg/ml (determined by Bradford assay)

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.