

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

DHH, 23-198aa, Human, His-tagged, Recombinant, E.coli

产品货号: TP01873

第三版

别名:desert hedgehog preproprotein, Drosophila, HHG-3, MGC35145, Desert hedgehog protein N-product

描述:DHH belongs to the hedgehog family. The C-terminal domain displays an autoproteolysis activity and a cholesterol transferase activity. The N- terminal domain is the active species in both local and longrange signaling, whereas the C-terminal domain has no signaling activity.

配方:Liquid in 20mM MES pH 5.5, 0.5mM DTT, 20% glycerol

分子量:22.0kDa (197aa)

序列:

MGSSHHHHHHSSGLVPRGSHMCGPGRGPVGRRRYARKQLVPLLYKQFVPGVPERTLGASGPAEGRVARGSERFRDLVP  
NYNPDIIIFKDEENSGADRLMTERCKERVNALAIVMNMWPGVRLRVTEGWDEDGHHAQDSLHYEGRALDITTSRDR  
NKYGLLARLAVEAGFDWVYYESRNHVHVSVKADNSLAVRAGG

纯度:> 95% by HPLC

浓度:1 mg/ml

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.