

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Dkk-1, 32-266aa, Human, His-tagged, Recombinant, Hi-5 cell

产品货号: TP01884

第三版

**别名:**Dickkopf-related protein 1, Dickkopf-1, hDkk-1, SK

**描述:**Dickkopf-related protein 1, also known as Dkk-1, is a member of the Dkk protein family including Dkk-1, 2, 3 and -4. It is a secreted protein with two cysteine rich regions and is involved in embryonic development through its inhibition of the WNT signaling pathway.

**配方:**Liquid. In PBS (pH 7.4) containing 30% glycerol, 2mM DTT, 2mM EDTA, 0.1mM PMSF.

**分子量:**27.8 kDa (253 aa), (On SDS-PAGE under denatured condition, apparent molecular weight of glycosylated rhDkk-1 protein will appear at approximately 40kDa)

**序列:**

ADPTLNSVLNSNAIKNLPPPLGGAAGHPGSAVSAAPGILYPGGNKYQTIDNYQPYPCAEDEECGTDEYCASPTRGGDAGV  
QICLACRKRRCMRHAMCCPGNYCKNGICVSSDQNHFRGEIETITESFGNDHSTLDGYSRRTLSSKMYHTKGQEGSV  
CLRSSDCASGLCCARHFWSKICKPVLKEGVCTKHRRKGSHGLEIFQRCYCGEGLSRIQKDHQASNSSRLHTCQRHSG  
RLVPRGSHHHHHH

**纯度:**> 95% by HPLC

**浓度:**0.25 mg/ml (determined by Bradford assay)

**内毒素:**<1.0 EU per 1 ug of protein (determined by LAL method)

**存储:** +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.