

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

DRS (Aspartyl-tRNA synthetase)

产品货号: TP01923

第三版

别名:

**描述:**Asparatyl tRNA synthetase is an alpha2 dimer that belongs to class IIb. Struntural analysis combined with mutagenesis and enzymology data on the yeast enzyme point to a tRNA binding process that starts by a recognition event between the tRNA anticodon loop and the synthetase anticodon binding module.

**配方:**1X PBS(pH7.2) 200mM NaCl 20% Glycerol

**分子量:**59kDa

**序列:**

**纯度:**> 95% by HPLC

**浓度:**0.50mg/ml

**内毒素:**<1.0 EU per 1 ug of protein (determined by LAL method)

**存储:** +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.