



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

DYNLRB2, 1-96aa, Human, His tag, E.coli

产品货号: TP01947

第三版

别名: Dynein light chain roadblock-type 2, DNCL2B, DNLC2B, MGC62033, ROBLD2

描述: DYNLRB2, also known as dynein light chain roadblock-type 2, acts as one of several non-catalytic accessory components of the cytoplasmic dynein complex that are thought to be involved in linking dynein to cargos and to adapter proteins that regulate dynein function.

配方: Liquid. 20mM Tris-HCl buffer (pH8.0) containing 20% glycerol 0.1M NaCl, 1mM DTT

分子量: 13.0 kDa (116aa) confirmed by MALDI-TOF

序列:

MGSSHHHHHSSGLVPRGSHMAEVEETLKRIQSHKGVIGTMVVNAEGPIRTTLDNSTTVQYAGLLHHLTMMKAKSTVRD
IDPQNDLTLRIRSKKHEIMVAPDKEYLLIVIQNPCE

纯度: > 95% by HPLC

浓度: 0.5 mg/ml (determined by Bradford assay)

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.