

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

EAF1, 182-268aa, Human, His tag, E.coli

产品货号: TP01952

第三版

别名: ELL associated factor 1, ELL-associated factor 1

描述: EAF1 belongs to the EAF family and is strongly expressed in heart, brain, placenta, lung, liver, skeletal muscle, kidney, pancreas, spleen, prostate, testis, small intestine and colon. This protein acts as a transcriptional transactivator of ELL and ELL2 elongation activities. Recombinant human EAF1 protein, fused to His-tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography techniques.

配方: Liquid. In 20mM Tris-HCl buffer (pH 8.0) containing 0.15M NaCl, 10% glycerol

分子量: 11.3 kDa (110aa) confirmed by MALDI-TOF

序列:

MGSSHHHHHHSSGLVPRGSHMGSDIIEQMSSSSGSSSSDSESSSGSDDSSSSGGEDNGPASPPQPSHQQPYNRPVANGT
SRPQGSNQLMNTLRNDLQLSESGSDSDD

纯度: > 95% by HPLC

浓度: 0.5 mg/ml (determined by BCA assay)

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.