
本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

EBAG9, 28-213aa, Human, His tagged, Recombinant, E.coli

产品货号: TP01953

第三版

别名: Estrogen receptor binding site associated, antigen 9, EB9, PDAF, RCAS1

描述: Estrogen receptor binding site associated antigen 9, also known as EBAG9, is a type III transmembrane protein that is predominantly located in the Golgi. This protein is a tumor-associated antigen that is expressed at high frequency in a variety of cancers, such as advanced breast and prostate cancers.

配方: Liquid in 20 mM Tris pH 8.0, 1 mM DTT, 2 mM EDTA, 10% glycerol

分子量: 23.4 kDa (207 aa) (Real molecular weight on SDS-PAGE will be shift up)

序列:

纯度: > 95% by HPLC

浓度: 1 mg/ml (determined by Bradford assay)

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.