



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

EDF1, 1-148aa, Human, His tag, E.coli

产品货号: TP01965

第三版

别名: Endothelial differentiation-related factor 1, EDF-1, MBF1

描述: Endothelial differentiation-related factor 1, also known as EDF1 may regulate endothelial cell differentiation. It has been postulated that the protein functions as a bridging molecule that interconnects regulatory proteins and the basal transcriptional machinery, thereby modulating the transcription of genes involved in endothelial differentiation.

配方: Liquid. In 20mM Tris-HCl buffer (pH 8.0) containing 1mM DTT, 20% glycerol, 0.1M NaCl

分子量: 17.4 kDa (156aa) confirmed by MALDI-TOF

序列:

MAESDWDTVTVLRKKGPAAQAKSKQAILAAQRGDEVETSKKWAAGQNKQHSITKNTAKLDRETEELHHDRVTLEV
GKVIQQGRQSKGLTQKDLATKINEKPQVIADYESGRAIPNNQVLGKIERAIGLKLKRGKDIGKPIEKGPRAKLEHHHHHH

纯度: > 95% by HPLC

浓度: 0.5 mg/ml (determined by Bradford assay)

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.