
本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human ERN1 / IRE1 Protein (aa 465-977)

产品货号: TP08589

第三版

别名:hIRE1p;IRE1;IRE1a;IRE1P

描述:Trypsin-3, also known as Trypsin III, brain trypsinogen, Serine protease 3 and PRSS3, is a secreted protein which belongs to the peptidase S1 family. Trypsin-3 / PRSS3 is expressed in pancreas and brain. It contains one peptidase S1 domain. Trypsin-3 / PRSS3 can degrade intrapancreatic trypsin inhibitors that protect against CP. Genetic variants that cause higher mesotrypsin activity might increase the risk for chronic pancreatitis (CP). A sustained imbalance of pancreatic proteases and their inhibitors seems to be important for the development of CP. The trypsin inhibitor-degrading activity qualified PRSS3 as a candidate for a novel CP susceptibility gene. Trypsin-3 / PRSS3 has been implicated as a putative tumor suppressor gene due to its loss of expression, which is correlated with promoter hypermethylation, in esophageal squamous cell carcinoma and gastric adenocarcinoma.

配方:PBS

分子量:58.3 kDa

序列:Pro 465-Leu 977

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.