



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Rat CLEC4D / CLECSF8 Protein (Fc tag)

产品货号: TP08592

第三版

别名:CLEC4D;Clecsf8;Mcl

描述:C-type lectin (CLEC) family is a type of carbohydrate-binding protein domain named lectin. C-type lectins are the most diverse and prevalent lectin family in immunity with its requirement for calcium for binding. Proteins including a C-type lectin domain have diverse range of functions including cell-cell adhesion, immune response to pathogens and apoptosis. There are at least 14 types of C-type lectins: type 1 to type 14. CLEC4D(CLECSF8) is a type 1 membrane glycoprotein belonging to the C-type lectin family, with restricted expression in the monocyte/macrophage lineage. It plays important roles in the function of macrophages.

配方:PBS

分子量:48.4 kDa

序列:Trp48-Lys218

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.