
本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human HIST3H2A / Histone H2A Protein

产品货号: TP08614

第三版

别名:MGC3165

描述: Histones are a complex family of highly conserved basic proteins responsible for packaging chromosomal DNA into nucleosomes. There are subtype diversities: H1, H2A, H2B and H3 or H4. It has become more and more evident that histone modifications are key players in the regulation of chromatin states and dynamics as well as in gene expression. Therefore, histone modifications and the enzymatic machineries that set them are crucial regulators that can control cellular proliferation, differentiation, plasticity, and malignancy processes. However, extracellular histones are a double-edged sword because they also damage host tissue and may cause death. Histones bound to platelets, induced calcium influx, and recruited plasma adhesion proteins such as fibrinogen to induce platelet aggregation. Histone cluster 3, H2a also known as histone H2A (HIST3H2A) is a member of histones. Covalent modification of histones is important in regulating chromatin dynamics and transcription. One example of such modification is ubiquitination, which mainly occurs on histones H2A and H2B. E3 ubiquitin ligase complex is specific for histone H2A (HIST3H2A). Reducing the expression of Ring2 results in a dramatic decrease in the level of ubiquitinated H2A in HeLa cells. DNA damage induces monoubiquitylation of histone H2A (HIST3H2A) in the vicinity of DNA lesions.

配方: PBS**分子量:** 14.2 kDa**序列:** Met 1-Lys 130**纯度:** > 95% by HPLC**浓度:****内毒素:** <1.0 EU per 1 ug of protein (determined by LAL method)**存储:** +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.