

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

TNFA, Donkey, His tag, E.coli

产品货号: TP08633

第三版

别名: Tumor necrosis factor alpha

描述: Tumor necrosis factor (TNF, cachexin, or cachectin; often called tumor necrosis factor alpha or TNF- α) is an adipokine and a cytokine. TNF is a member of the TNF superfamily, which consists of various transmembrane proteins with a homologous TNF domain.

配方: 10mM PBS, pH7.4, containing 1mM DTT, 5% trehalose.

分子量: 27.0kDa

序列: Met1~Leu234

纯度: > 95% by HPLC

浓度: 200ug/ml

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.