

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

LTF,Donkey,His tag ,E.coli

产品货号: TP08634

第三版

别名:Lactoferrin

描述:Lactoferrin, also known as lactotransferrin, is a multifunctional protein of the transferrin family. Lactoferrin is a globular glycoprotein with a molecular mass of about 80 kDa that is widely represented in various secretory fluids, such as milk, saliva, tears, and nasal secretions. Lactoferrin is also present in secondary granules of PMNs and is secreted by some acinar cells. Lactoferrin can be purified from milk or produced recombinantly. Human colostrum has the highest concentration, followed by human milk, then cow milk.

配方:10mM PBS, pH7.4, containing 1mM DTT, 5% trehalose.

分子量:9.0kDa

序列:Met1~Phe23

纯度:> 95% by HPLC

浓度:200ug/ml

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.