

**14-3-3 ? Protein (63A1) Anti-Human Mouse IgG MoAb**

**Monoclonal Antibody**

**Cat.NO.: MA01017**

---

3th Edition

**Immunogen:** Synthetic peptide in portion of C-terminus of Human 14-3-3 ? (QQDDDGEGNN)

**Isotype:** Mouse-Mouse hybridoma (Supernatant) IgG1

**Clone:** 63A1

**Description:**

**Form:** Lyophilized product from 1% BSA in PBS containing 0.05% NaN<sub>3</sub>

**Concentration:** 100 ug/ml

**Storage:** Can be stored at +4°C. For long term storage, aliquot and store at -20°C. Avoid repeated freezing and thawing cycles.

**Usage:** The antibody has been tested by ELISA, Western blot analysis and ICC/IF to assure specificity and reactivity. Since application varies, however, each investigation should be titrated by the reagent to obtain optimal results.

**Application:** ELISA, WB, ICC/IF