

Amyloid Beta (35-40) (1A10) A? Anti-Human Mouse IgG MoAb

Monoclonal Antibody

Cat.NO.: MA01078

3th Edition

Immunogen: Synthetic peptide of human Amyloid Beta A? (35-40) (MVGGVV).

Isotype: IgG1

Clone: 1A10

Description:

Form: Lyophilized product in PBS containing 1 % BSA and 0.05 % NaN₃.

Concentration: 1.0 mL deionized water will be added to the product, then its concentration comes to 50 µg/mL.

Storage: Can be stored at +4°C. For long term storage, aliquot and store at -20°C. Avoid repeated freezing and thawing cycles.

Usage: The antibody has been tested by ELISA, Western blot analysis and ICC/IF to assure specificity and reactivity. Since application varies, however, each investigation should be titrated by the reagent to obtain optimal results.

Application: ELISA, WB, ICC/IF