

CD20 (N) (4-6H:2C) Anti-Human Mouse IgG MoAb

Monoclonal Antibody

Cat.NO.: MA01151

3th Edition

Immunogen: Synthetic peptide of the N-terminal part of human CD20 (CMQSGPKPLFRRMSS).

Isotype: IgG3

Clone: 4-6H:2C

Description:

Form: Lyophilized product from 1.0 mL of PBS containing 1 % BSA and 0.05 % NaN₃

Concentration: Please see datasheet.

Storage: Can be stored at +4°C. For long term storage, aliquot and store at -20°C. Avoid repeated freezing and thawing cycles.

Usage: The antibody has been tested by ELISA, Western blot analysis and ICC/IF to assure specificity and reactivity. Since application varies, however, each investigation should be titrated by the reagent to obtain optimal results.

Application: ELISA, WB, ICC/IF