

Anti-CD252/OX40L (TAG34) Mouse IgG MoAb

Monoclonal Antibody

Cat.NO.: MA01156

3th Edition

Immunogen:Lysate of HTLV-1 persistently-infected cell line (MT-2)

Isotype:IgG1

Clone:TAG34

Description:

Form:Lyophilized product from PBS containing 1 % BSA and 0.05 % NaN₃

Concentration:Please see datasheet.

Storage: Can be stored at +4°C. For long term storage, aliquot and store at -20°C. Avoid repeated freezing and thawing cycles.

Usage:The antibody has been tested by ELISA, Western blot analysis and ICC/IF to assure specificity and reactivity. Since application varies, however, each investigation should be titrated by the reagent to obtain optimal results.

Application:ELISA, WB, ICC/IF