

GARS Monoclonal Antibody (Clone AT4E10)

Monoclonal Antibody

Cat.NO.: MA01272

3th Edition

Immunogen: Recombinant human GARS (43-289aa) purified from E. coli

Isotype: IgG1,k

Clone: Anti-human GARS mAb, clone AT4E10, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human GARS protein.

Description:

Form: Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

Concentration: 1 mg/mL

Storage: Can be stored at +4°C. For long term storage, aliquot and store at -20°C. Avoid repeated freezing and thawing cycles.

Usage: The antibody has been tested by ELISA, Western blot analysis and ICC/IF to assure specificity and reactivity. Since application varies, however, each investigation should be titrated by the reagent to obtain optimal results.

Application: ELISA, WB, ICC/IF