

H1N1/HA1 Anti-Human Antibody (Clone AT1G7)

Monoclonal Antibody

Cat.NO.: MA01302

3th Edition

Immunogen: Recombinant human H1N1/HA1 (18-344aa) purified from Baculovirus

Isotype: Mouse IgG1 heavy chain and ? light chain

Clone: 1G7

Description:

Form: Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

Concentration: 1mg/ml

Storage: Can be stored at +4°C. For long term storage, aliquot and store at -20°C. Avoid repeated freezing and thawing cycles.

Usage: The antibody has been tested by ELISA, Western blot analysis and ICC/IF to assure specificity and reactivity. Since application varies, however, each investigation should be titrated by the reagent to obtain optimal results.

Application: ELISA, WB, ICC/IF