

**H5N1/HA1 Anti-Human Antibody (Clone AT2B7)**

**Monoclonal Antibody**

**Cat.NO.: MA01304**

---

3th Edition

**Immunogen:** Recombinant human H5N1/HA1 (17-338aa) purified from Baculovirus

**Isotype:** Mouse IgG1 heavy chain and kappa light chain

**Clone:** 2B7

**Description:**

**Form:** Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

**Concentration:** 1mg/ml

**Storage:** Can be stored at +4°C. For long term storage, aliquot and store at -20°C. Avoid repeated freezing and thawing cycles.

**Usage:** The antibody has been tested by ELISA, Western blot analysis and ICC/IF to assure specificity and reactivity. Since application varies, however, each investigation should be titrated by the reagent to obtain optimal results.

**Application:** ELISA, WB, ICC/IF