

**Hexokinase 3 Anti-Human Antibody (Clone 4F5)**

**Monoclonal Antibody**

**Cat.NO.: MA01308**

---

3th Edition

**Immunogen:** Recombinant human Hexokinase 3 (1-923aa) purified from E. coli

**Isotype:** Mouse IgG1 heavy chain and kappa light chain

**Clone:** Anti-human Hexokinase 3 mAb, clone 4F5, is derived from hybridization of mouse SP2/0 myeloma cells with spleen cells from BALB/c mice

**Description:**

**Form:** Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

**Concentration:** 1 mg/ml

**Storage:** Can be stored at +4°C. For long term storage, aliquot and store at -20°C. Avoid repeated freezing and thawing cycles.

**Usage:** The antibody has been tested by ELISA, Western blot analysis and ICC/IF to assure specificity and reactivity. Since application varies, however, each investigation should be titrated by the reagent to obtain optimal results.

**Application:** ELISA, WB, ICC/IF