



**Tenascin-C (EGF Like Domain) (4F10TT) Anti-Human Mouse IgG**

**Monoclonal Antibody**

**Cat.NO.: MA01329**

---

3th Edition

**Immunogen:**Purified human glioma tenascin-C

**Isotype:**IgG1

**Clone:**4F10TT

**Description:**

**Form:**Lyophilized product from 1% BSA in PBS containing 0.05%NaN3

**Concentration:**Please see datasheet.

**Storage:** Can be stored at +4°C. For long term storage, aliquot and store at -20°C. Avoid repeated freezing and thawing cycles.

**Usage:**The antibody has been tested by ELISA, Western blot analysis and ICC/IF to assure specificity and reactivity. Since application varies, however, each investigation should be titrated by the reagent to obtain optimal results.

**Application:**ELISA, WB, ICC/IF