

Anti-Human/Mouse/Rat PI 3 kinase p85 alpha Monoclonal Antibody

Monoclonal Antibody

Cat.NO.: MA02087

3th Edition

Immunogen:Fusion protein of PI3K p85(alpha)

Isotype:IgG2a

Clone:645

Description:

Form:

Concentration:None

Storage: Can be stored at +4°C. For long term storage, aliquot and store at -20°C. Avoid repeated freezing and thawing cycles.

Usage:The antibody has been tested by ELISA, Western blot analysis and ICC/IF to assure specificity and reactivity. Since application varies, however, each investigation should be titrated by the reagent to obtain optimal results.

Application:ELISA, WB, ICC/IF