



**Anti-Human/Mouse/Rat PI 3 kinase p85 alpha Monoclonal Antibody**

**Monoclonal Antibody**

**Cat.NO.: MA02087**

---

3th Edition

**Immunogen:**Fusion protein of PI3K p85(alpha)

**Isotype:**IgG2a

**Clone:**645

**Description:**

**Form:**

**Concentration:**None

**Storage:** Can be stored at +4°C. For long term storage, aliquot and store at -20°C. Avoid repeated freezing and thawing cycles.

**Usage:**The antibody has been tested by ELISA, Western blot analysis and ICC/IF to assure specificity and reactivity. Since application varies, however, each investigation should be titrated by the reagent to obtain optimal results.

**Application:**ELISA, WB, ICC/IF