



Phosphoserine Monoclonal Antibody

Monoclonal Antibody

Cat.NO.: MA02606

3th Edition

Immunogen:Purified Protein

Isotype:IgG

Clone:7A3

Description:

Form:

Concentration:1mg/mL

Storage: Can be stored at +4°C. For long term storage, aliquot and store at -20°C. Avoid repeated freezing and thawing cycles.

Usage:The antibody has been tested by ELISA, Western blot analysis and ICC/IF to assure specificity and reactivity. Since application varies, however, each investigation should be titrated by the reagent to obtain optimal results.

Application:ELISA, WB, ICC/IF