

Anti-Human/Mouse/Rat Phospho-ERK 1/2 (Tyr222/205) Monoclonal Antibody

Monoclonal Antibody

Cat.NO.: MA02627

3th Edition

Immunogen: Synthetic Peptide of Phospho-ERK 1/2 (Tyr222/205)

Isotype: IgG

Clone: 4D4

Description:

Form:

Concentration: 1mg/mL

Storage: Can be stored at +4°C. For long term storage, aliquot and store at -20°C. Avoid repeated freezing and thawing cycles.

Usage: The antibody has been tested by ELISA, Western blot analysis and ICC/IF to assure specificity and reactivity. Since application varies, however, each investigation should be titrated by the reagent to obtain optimal results.

Application: ELISA, WB, ICC/IF