



Anti-Human/Mouse/Rat MAP2K2 Monoclonal Antibody

Monoclonal Antibody

Cat.NO.: MA02656

3th Edition

Immunogen: Synthesized peptide derived from human MEK-2 around the non-phosphorylation site of T394.

Isotype: IgG

Clone:

Description:

Form:

Concentration: 1mg/mL

Storage: Can be stored at +4°C. For long term storage, aliquot and store at -20°C. Avoid repeated freezing and thawing cycles.

Usage: The antibody has been tested by ELISA, Western blot analysis and ICC/IF to assure specificity and reactivity. Since application varies, however, each investigation should be titrated by the reagent to obtain optimal results.

Application: ELISA, WB, ICC/IF