

### **14-3-3 ? Protein Anti-Human Rabbit IgG Affinity Purify**

**Polyclonal Antibody**

**Cat.NO.: PA01024**

---

3th Edition

**Description:** The 14-3-3 proteins are a family of conserved regulatory molecules expressed in all eukaryotic cells. A striking feature of the 14-3-3 proteins is their ability to bind a multitude of functionally diverse signaling proteins, including kinases, phosphatases, and transmembrane receptors.

**Antigen:** Synthetic peptide for N-terminal of Human 14-3-3 ?.

**Form:** Lyophilized product from 1% BSA in PBS containing 0.05% NaN<sub>3</sub>.

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:** Human 14-3-3 ? specific, Non-crossreact to Human 14-3-3 ?, ?, ?, ?, ?.