

Angiotensinogen (405) Anti-Mouse/Rat Rabbit IgG Affinity Purify

Polyclonal Antibody

Cat.NO.: PA01050

3th Edition

Description:Angiotensinogen is the precursor of angiotensin and is cleaved into angiotensin I and II in the renin-angiotensin system, and it has long been reported to play an important role in controlling blood pressure.

Antigen:Synthetic peptide of the C terminal part of Mouse Angiotensinogen (IGDTNPRVGEVLNSIL).

Form:Lyophilized product from 1% BSA in PBS containing 0.05% NaN₃.

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Reacts with mouse and rat Angiotensinogen.