

ApoE (A299) Anti-Human Rabbit IgG Affinity Purify

Polyclonal Antibody

Cat.NO.: PA01052

3th Edition

Description: Apolipoprotein-E (Apo-E) is a structural component of very low-density lipoprotein (VLDL) synthesized by the liver and intestinally synthesized chylomicrons.

Antigen: Synthetic peptide for c-terminal of Human ApoE.

Form: Lyophilized product from 1% BSA in PBS containing 0.05% NaN₃.

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Can be react to ApoE2, ApoE3 and ApoE4.