

**Tie (N1125) Anti-Human Rabbit IgG Affinity Purify**

**Polyclonal Antibody**

**Cat.NO.: PA01202**

---

3th Edition

**Description:** This antibody can be stained in formalin fixed paraffin embedded tissues after microwave pretreatment by several Immunohistochemical techniques such as Avidin Biotin Complex (ABC) method.

**Antigen:** Part of synthetic peptide of Human Tie-1

**Form:** Lyophilized product from 1% BSA in PBS containing 0.05% NaN<sub>3</sub>

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:** Confirmed by Tie/BaF3 Transfectant and HEL