

TNF- α Anti-Rat Rabbit IgG Affinity Purify

Polyclonal Antibody

Cat.NO.: PA01205

3th Edition

Description: TNF α (Tumor Necrosis Factor alpha) is a 17.5 kDa, 157 amino acid protein that is a potent lymphoid factor, which exerts cytotoxic effects on a wide range of tumor cells and other target cells.

Antigen: Synthetic peptide for C-Terminal of rat TNF α

Form: Lyophilized product from 1% BSA in PBS containing 0.05%NaN₃

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1 μ g/ml.

Specificity: