

**VEGF (V-N) Anti-Human Rabbit IgG Affinity Purify**

**Polyclonal Antibody**

**Cat.NO.: PA01215**

---

3th Edition

**Description:** Vascular Endothelial Cell Growth Factor (VEGF) is a homodimeric protein initially purified from media conditioned by normal bovine pituitary folliculo-stellate cells and secreted by a variety of vascularized tissues

**Antigen:** Synthetic peptide for N-terminal of VEGF

**Form:** Lyophilized product from 1% BSA in PBS containing 0.05% NaN<sub>3</sub>

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1 µg/ml.

**Specificity:** Human specific. Cross reactivity with Mice and Rat are not tested. React with VEGFs both of reduced form and of non-reduced form.