

VEGF (V-N) Anti-Rat Rabbit IgG Affinity Purify

Polyclonal Antibody

Cat.NO.: PA01217

3th Edition

Description: Vascular Endothelial Cell Growth Factor (VEGF) is a homodimeric protein initially purified from media conditioned by normal bovine pituitary folliculo-stellate cells and secreted by a variety of vascularized tissues.

Antigen: Synthetic peptide for N-terminal of [Rat VEGF

Form: Lyophilized product from 1% BSA in PBS containing 0.05%NaN₃

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Rat specific. Cross react with Human VEGF (165). React with rat VEGFs both of reduced form and of non-reduced form.