

XIAP Anti-Human Rabbit IgG Affinity Purify

Polyclonal Antibody

Cat.NO.: PA01222

3th Edition

Description: IAP (inhibitor of apoptosis protein) is an evolutionarily-conserved inhibitory factor of apoptosis. IAP directly binds to caspase which plays the central role in operation of apoptosis and inhibits it.

Antigen: Synthetic peptide of a part of Human XIAP (aa244-263: VSSDRNFPNSTNLPRNPSMA)

Form: Lyophilized product from 1 % BSA in PBS containing 0.05% NaN₃

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Cross reacts with mouse.