

**GRO/CINC-2?, 2? (N) Anti-Rat Rabbit IgG Affinity Purify**

**Polyclonal Antibody**

**Cat.NO.: PA01232**

---

3th Edition

**Description:** Cytokine-induced neutrophil chemo attractant 1 (CINC-1) was originally purified from media conditioned by IL-1  $\alpha$  stimulated rat kidney epithelioid cells (NRK-52E).

**Antigen:** Synthetic peptide for N-Terminal common sequences of rat GRO/CINC-2? and rat GRO/CINC-2? conjugating bovine thyroglobulin.

**Form:** Lyophilized product from 1% BSA in PBS containing 0.05% NaN<sub>3</sub>

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:** React with Rat GRO/CINC-2 ? and –2 ? . Non-cross react with Rat GRO/CINC-1 and Rat GRO/CINC-3.