

**hnRNP B1 Anti-Human Rabbit IgG Affinity Purify**

**Polyclonal Antibody**

**Cat.NO.: PA01252**

---

3th Edition

**Description:** hnRNP proteins are a family at least 30 polypeptides, play a role in the splicing of pre-mRNA into the mature product in association with small nuclear ribonucleoprotein (snRNP) and splicing factors.

**Antigen:** Synthetic peptide of part of N terminal of human hnRNP (KTLETVPLERKKREK)

**Form:** Lyophilized product from 1% BSA in PBS containing 0.05% NaN<sub>3</sub>

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:** Confirmed by western blotting (37 kD).