

**KC Anti-Mouse Rabbit IgG**

**Polyclonal Antibody**

**Cat.NO.: PA01277**

---

3th Edition

**Description:** Mouse KC (Keratinocyto-derived Cytokine) is a counterpart of human GRO protein. It is reported that KC is produced in such as activated monocyte and fibroblast.

**Antigen:** Recombinant Mouse KC

**Form:** Lyophilized product from 1% BSA in PBS containing 0.05% NaN<sub>3</sub>

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**