

Leptin (C) Anti-Mouse Rabbit IgG Affinity Purify

Polyclonal Antibody

Cat.NO.: PA01283

3th Edition

Description:Leptin is a 16 KDa secreted protein produced by the ob gene. Adipose tissue produces Leptin and releases it into the bloodstream. As fat deposits grow, blood Leptin levels tends to increase.

Antigen:Synthetic peptide in portion of C-terminus of Mouse Leptin

Form:Lyophilized product from 1% BSA in PBS containing 0.05% NaN₃

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Cross-reacts with Rat