

LGP2 Anti-Human Rabbit IgG Affinity Purify

Polyclonal Antibody

Cat.NO.: PA01285

3th Edition

Description: A full-length cDNA encoding mouse LGP2 was cloned by computational analysis, RT-PCR, and RACE of liver and mammary tissue total RNA. And human LGP2 was also cloned.

Antigen: Synthetic peptide of part of c terminal of human LGP2 (IQAKKWSRVPFSV)

Form: Lyophilized product from 1% BSA in PBS containing 0.05% NaN₃

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Reacts with human and mouse LGP2