

**Anti-Mouse IgG (H+L) Goat IgG Affinity Purify**

**Polyclonal Antibody**

**Cat.NO.: PA01312**

---

3th Edition

**Description:** Anti-Mouse IgG (H+L) Goat IgG Affinity Purify, 1 mg. Lyophilized product from 1% BSA in PBS containing 0.05% NaN<sub>3</sub>, purified with antigen.

**Antigen:** Mouse IgG (H+L)

**Form:** Lyophilized product from 1% BSA in PBS containing 0.05% NaN<sub>3</sub>.

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:** The antibody is absorbed by rat IgG, goat IgG, rabbit IgG and human IgG. Cross-reactivity to goat IgG, rabbit IgG and human IgG is less than 0.1 %. Cross-reactivity to rat IgG is about 0.5 %.